



Glass Ampoules




Form B
Volumes: 1 ml to 50 ml
Break systems: OPC, CBR, SCO (up to 30 ml)




Form C
Volumes: 1 ml to 50 ml
Break systems: OPC, CBR, SCO (up to 30 ml)



Form D
Volumes: 1 ml to 30 ml
Break systems: OPC, CBR, SCO



Double tip
Volumes: 1 ml to 20 ml
Break systems: CBR, DRB



Marzocchi and Special Designs
Volumes: 1.5 ml to 20 ml
Break systems: OPC, CBR



Fine tip
Volumes: 1 ml to 15 ml
Break systems: CBR, SCO

Glass Vials



Injection vials for liquid formulations

- 1-100 ml
- StandardLine, ClearLine, TopLine
- FIOLAX® or DURAN® glass



Injection vials for dry powder formulations

- 1-100 ml
- StandardLine, ClearLine, TopLine
- FIOLAX® or DURAN® glass



Screw neck vials

- 1-100 ml
- StandardLine, ClearLine, TopLine
- FIOLAX® or DURAN® glass



Double chamber vials

- 1 ml, 2 ml, 4 ml, 8 ml, 10 ml and custom sizes
- StandardLine and additional specifications
- FIOLAX® or DURAN® glass



Tablet vials

- 1-100 ml
- StandardLine and additional specifications
- ILLAX® glass



Perfume vials

- Custom sizes up to 50 ml
- Silkscreen printing with heavy-metal-free precious or standard enamels



- 1-100 ml
- ClearLine, TopLine
- Barrier coating – resistant to delamination




- 1-10 ml
- ClearLine, TopLine
- Hydrophobic coating and improved bottom geometry



- 1-10 ml
- ClearLine, TopLine
- Hydrophobic coating





Sterile syringe sets – forma 3^s®



The delivery in nest form enables easy handling and filling of the syringes.

forma 3^s® product portfolio

- 0.5 ml to 3 ml
- Staked needle, luer cone, luer lock, luer lock V-OVS^{®1}
- Cut flange or round flange
- Printing or graduation available
- Other formats are available on demand

					
 Staked Needle 0.5ml, 1ml long, 1ml short, 22G – 29G	✓	✓			
 Luer Cone 1ml long, 1-3ml					✓
 Luer Lock 1ml long, 1-3ml			✓	✓	

Glass Cartridges



Pen cartridges

- Recommended designs from 1 to 5 ml
- ClearLine, TopLine



By-pass cartridges

- Customized sizes
- ClearLine, TopLine



Customized cartridges

- 1 to 10 ml
- ClearLine, TopLine



Dental cartridges

- 1.8 ml, 2.2 ml
- StandardLine